

Adding a workflow to BIAFLOWS




neubiasacademy.org



Defragmentation:
bringing BioImage Analysts to the cloud!



Overview




- Prerequisites:
 - A github account
 - A dockhub account
 - configure BIAFLOWS to scan your github repositories
- Create a new repository on github
 - Give the new repository the right to access your dockerhub
- Modify your workflow to
 - Run on the images in a folder
 - Accept input-parameters
 - Produce results in the expected form
- Add 4 files to the github repository
 - Dockerfile
 - Defines the computational environment
 - wrapper.py
 - Download images, run workflow, create/upload annotations, create/upload metrics
 - macro.ijm
 - the image analysis workflow
 - descriptor.json
 - Definition of the parameters of the workflow
- Make a release
 - A release triggers a git-hub action that
 - creates a docker-image from the repository
 - Pushes the docker-image to dockerhub

Cloud components



Poll for new docker images

Link from workflow to the repository



On release: create a docker image
Push it to dockerhub

Download and run a docker image



Example

ImageJ – 3D Spot Detection

- Basic idea
 - Apply Top-Hat-Transform
 - Find maxima
- Parameters:
 - Radius median in xy
 - Radius median in z
 - Radius top-hat in xy
 - Radius top-gat in z
 - Dynamic (ext. maxima)
 - 3D connectivity 6 or 26

Clear first and last slice
of the stack

Median 3D

White Top Hat

Ext. Maxima

Connected Components

Analyze Regions 3D

Example – The original macro

```
noise_filter_radius_xy = 1;  
noise_filter_radius_z = 0.5;  
top_hat_radius_xy = 7.5;  
top_hat_radius_z = 3.5;  
dynamic = 75;  
connectivity = 6;
```

```
getDimensions(width, height, channels, slices, frames);
```

```
Stack.setSlice(1);
```

```
run("Select All");
```

```
run("Clear", "slice");
```

```
Stack.setSlice(slices);
```

```
run("Clear", "slice");
```

```
run("Select None");
```


```
run("Median 3D...", "x="+noise_filter_radius_xy+" y="+noise_filter_radius_xy+" z="+noise_filter_radius_z);
```

```
run("Morphological Filters (3D)", "operation=[White Top Hat] element=Cube x-radius="+top_hat_radius_xy+" y-radius="+top_hat_radius_xy+" z-radius="+top_hat_radius_z);
```

```
run("Extended Min & Max 3D", "operation=[Extended Maxima] dynamic="+dynamic+" connectivity="+connectivity);
```

```
run("Connected Components Labeling", "connectivity="+connectivity+" type=[16 bits]");
```

```
run("Analyze Regions 3D", "centroid surface_area_method=[Crofton (13 dirs.)] euler_connectivity=6");
```



dockerhub


- Docker hub hosts and serves built docker-images
- You need a dockerhub account with the same name as your github account
 - dockerhub allows only letters and digits
 - '-' on github will be removed for dockerhub
 - everything will be made lower-case for dockerhub

volker-baecker volkerbaecker

- Organizations can also be used on both sites

Neubias-WG5 neubiaswg5


- Create a dockerhub account that corresponds to your github account!




Configure BIAFLOWS to scan your repositories

- In BIAFLOWS
 - Go to Workflows

- Fill in
 - Your github user-name
 - Your dockerhub user-name
 - The prefix of workflows to import



- Add a new trusted source





Source code provider	<input type="text" value="github"/>
Source code provider username	<input type="text" value="volker-baecker"/>
Source code provider token (optional)	<input type="text"/>
Environment provider	<input type="text" value="docker"/>
Environment provider username	<input type="text" value="volkerbaecker"/>
Prefix	<input type="text" value="W_"/>

Create a new repository on github

- Look for the most similar NEUBIAS-WG5 workflow repository

<https://github.com/orgs/Neubias-WG5/repositories>





- Press the “Use this template” button
- Add a name for the new repository, starting with W_
- Modify the Readme

Create a new repository from W_SpotDetection3D-IJ

The new repository will start with the same files and folders as [Neubias-WG5/W_SpotDetection3D-IJ](#).

Owner *



volker-baecker ▾

Repository name *

W_SpotDetection3D-IJ-TopHat



Great repository names are short and memorable. Need inspiration? How about **expert-palm-tree**?

Description (optional)



Public

Anyone on the internet can see this repository. You choose who can commit.



Private

You choose who can see and commit to this repository.








Include all branches

Copy all branches from Neubias-WG5/W_SpotDetection3D-IJ and not just master.



You are creating a public repository in your personal account.

Create repository from template

 .github/workflows	Initial commit	4 minutes ago
 Dockerfile	Initial commit	4 minutes ago
 IJSpotDetection3D.ijm	Initial commit	4 minutes ago
 README.MD	Update README.MD	3 minutes ago
 descriptor.json	Initial commit	4 minutes ago
 wrapper.py	Initial commit	4 minutes ago

README.MD



W_SpotDetection3D-IJ-TopHat

This workflow detects spots in a 3D image by filtering the image by a Top-Hat-Filter (user defined radius) and detecting local intensity maxima with a given dynamic (user defined threshold).

Give the new repository the right to access your dockerhub

- Log in to dockerhub
- Click on your username
 - In the upper right corner
 - Go to “Account Settings”
- Select Security
- Press “New Access Token”
- Enter a name for the token
- Generate the token
- **Copy the token**

New Access Token

A personal access token is similar to a password except you can have many tokens and revoke access to each one at any time. [Learn more](#)

Access Token Description *

github-token

Access permissions


Read, Write, Delete



Read, Write, Delete tokens allow you to manage your repositories.

Give the new repository the right to access your dockerhub

- Go back to your github-repository
- In “Settings”
 - Go to “Secrets>Actions”
- Create a secret with the name DOCKERHUB_TOKEN and paste the copied token as its value
- Create a secret with the name DOCKERHUB_USERNAME and your dockerhub username as value
- Instead of a doing it for each repository you can do it for an organization



Adapt your workflow

- Replace the macro in the git-repository with your macro
- Modify it to produce the expected result
 - Mask of centroids
- Modify the macro to accept parameters from the command-line
- Modify it to work on all images in the input folder and write the results to the output folder

Create the expected result

- What is the expected result?
- You find the information in the BIAFLOWS documentation, see:
- In our case
 - A 16-bit image
 - With a pixel of 65535 for each centroid of a spot
 - 0 for the background

https://neubias-wg5.github.io/problem_class_ground_truth.html#steps-section

Create the expected result

- Our macro just measures the centroids of the spots and reports them in a table
- From the table we create the output image, by setting a pixel to 65535 for each centroid in the table

```
columnNames = split(Table.headings, "\t");
X = Table.getColumn(columnNames[1]);
Y = Table.getColumn(columnNames[2]);
Z = Table.getColumn(columnNames[3]);
Table.sort(columnNames[3]);
run("Select All");
run("Clear", "stack");
run("Select None");
for (c = 0; c < X.length; c++) {
    x = X[c];
    y = Y[c];
    z = Z[c];
    Stack.setSlice(z);
    setPixel(x, y, 65535);
}
```


Read parameters from the command-line

- Our macro has the parameters
 - noise_filter_radius_xy
 - noise_filter_radius_z
 - top_hat_radius_xy
 - top_hat_radius_z
 - dynamic
 - Connectivity
- In addition we need parameters for the input and output folders
 - inputDir
 - outputDir

```
// Path to input image and output mask
inputDir = "/media/baecker/DONNEES/mri/2022/neubias/data";
outputDir = "/media/baecker/DONNEES/mri/2022/neubias/out";

// Functional parameters
noise_filter_radius_xy = 1;
noise_filter_radius_z = 0.5;
top_hat_radius_xy = 7.5;
top_hat_radius_z = 3.5;
dynamic = 75;
connectivity = 6;

arg = getArgument();
parts = split(arg, ",");

for(i=0; i<parts.length; i++) {
    nameAndValue = split(parts[i], "=");
    if (indexOf(nameAndValue[0], "input") > -1) inputDir = nameAndValue[1];
    if (indexOf(nameAndValue[0], "output") > -1) outputDir = nameAndValue[1];
    if (indexOf(nameAndValue[0], "filter_xy") > -1) noise_filter_radius_xy = nameAndValue[1];
    if (indexOf(nameAndValue[0], "filter_z") > -1) noise_filter_radius_z = nameAndValue[1];
    if (indexOf(nameAndValue[0], "top_hat_xy") > -1) top_hat_radius_xy = nameAndValue[1];
    if (indexOf(nameAndValue[0], "top_hat_z") > -1) top_hat_radius_z = nameAndValue[1];
    if (indexOf(nameAndValue[0], "dynamic") > -1) dynamic = nameAndValue[1];
    if (indexOf(nameAndValue[0], "connectivity") > -1) connectivity = nameAndValue[1];
}
```

Run on all images in a folder

- We got the input and output folders as parameters
- Apply workflow to each image in the input folder and write the result to the output folder
- Cleanup after each iteration
- Run in batch-mode
- Quit the java-virtual machine process at the end

```
setBatchMode(true);

images = getFileList(inputDir);

for(i=0; i<images.length; i++) {
    image = images[i];
    if (!endsWith(image, ".tif")) continue;
    open(inputDir + "/" + image);

    ...

    // Export results
    save(outputDir + "/" + image);

    // Cleanup
    run("Close All");
    close(title + "-morpho");
}
run("Quit");
```

Setting up the environment

- We have to install the ImageJ-plugin MorphoLibJ
- Use a more recent version of FIJI

We change the Dockerfile accordingly

```
# Install Fiji plugins
```

```
RUN cd /fiji/plugins && wget -O MorphoLibJ_-1.5.1.jar
```

```
https://github.com/ijpb/MorphoLibJ/releases/download/v1.5.1/MorphoLibJ\_-1.5.1.jar
```

```
# Install Fiji.
```

```
RUN wget --no-check-certificate https://downloads.imagej.net/fiji/archive/20221013-0917/fiji-linux64.zip
```

```
RUN unzip fiji-linux64.zip
```

```
...
```

```
# Clean up
```

```
RUN rm fiji-linux64.zip
```

Tell the wrapper how to run your workflow

- The wrapper downloads images, runs your workflow and uploads results
- The only thing we need to change is the passing of the parameters from BIAFLOWS GUI to your workflow

```
command = "/usr/bin/xvfb-run java -Xmx6000m -cp /fiji/jars/ij.jar ij.ImageJ --headless --console " \
    "-macro macro.ijm \"input={}, output={}, filter_xy={}, filter_z={}, top_hat_xy={}, \
top_hat_z={}, dynamic={}, connectivity={}\" \" \" \
    .format(in_path, out_path, nj.parameters.filter_xy, nj.parameters.filter_z, \
nj.parameters.top_hat_xy, nj.parameters.top_hat_z, nj.parameters.dynamic, \
nj.parameters.connectivity)
```


Adapt the json descriptor

- Fields to be changed
 - Name
 - Image (on dockerhub)
 - Description
 - Command-line (parameters)
 - Inputs (parameters)
 - after cytomine_id_software
 - Set GUI text and default values for the parameters

```
"name": "SpotDetection3D-IJ-TopHat",  
  
"image": "volkerbaecker/w_spotdetection3d-ij-tophat",  
  
"description": "3D spot detection using a top hat filter and the detection of maxima.",  
  
"command-line": "python wrapper.py CYTOMINE_HOST CYTOMINE_PUBLIC_KEY  
CYTOMINE_PRIVATE_KEY CYTOMINE_ID_PROJECT CYTOMINE_ID_SOFTWARE  
FILTER_XY FILTER_Z TOP_HAT_XY TOP_HAT_Z, DYNAMIC, CONNECTIVITY",  
  
{  
    "id": "filter_xy",  
    "value-key": "@ID",  
    "command-line-flag": "--@id",  
    "name": "Radius xy",  
    "description": "Radius for the Median-filter in the x-y-plane",  
    "type": "Number",  
    "default-value": 1,  
    "optional": true  
}, ...
```

Create a release on github


- Click on tags
- Create a new release
- Use a version number of the form
 $v\langle x \rangle.\langle y \rangle.\langle z \rangle$



The screenshot shows the GitHub 'Create a new release' page. At the top, there are two tabs: 'Releases' (active) and 'Tags'. Below the tabs, there are two dropdown menus: 'v0.1' and 'Target: master'. A message below these says 'Excellent! This tag will be created from the target when you publish this release.' Below this is a text input field labeled 'Release title'. Underneath the input field are two tabs: 'Write' (active) and 'Preview'. Below the tabs is a rich text editor with various formatting icons (bold, italic, link, etc.). The text '3D spot detection with top-hat filter and finding maxima' is entered in the editor.


Check if the docker image is built correctly

- In github click on “Actions”
- Look for the version of the release you just created
 - Note: Building the image can take some time



Add the workflow to a problem class in BIAFLOWS

- Login to BIAFLOWS
- Select a problem
 - Problem: SPOT-COUNTING-3D
- Click on configuration
- Click on workflows
- Search your workflow and enable it



The screenshot shows a web interface for managing workflows. At the top, there is a search bar containing the text 'SpotDetection3D'. Below the search bar is a table with four columns: 'Name', 'Version', 'Runnable', and 'Status'. The table lists four workflows. The first three are 'SpotDetection3D-Icy (v2.1)', 'SpotDetection3D-Hessian-IJ (v1.2)', and 'SpotDetection3D-IJ (1.7)', all with 'Last release' versions, 'Runnable' status, and 'Enabled' status. The fourth workflow is 'SpotDetection3D-Icy (v2.0)', which is 'Deprecated', 'Runnable', and 'Disabled'. At the bottom left, there is a '10 per page' dropdown menu. At the bottom right, there is a pagination control showing '1' and navigation arrows.

Name	Version	Runnable	Status
SpotDetection3D-Icy (v2.1)	Last release	✓ Yes	Enabled
SpotDetection3D-Hessian-IJ (v1.2)	Last release	✓ Yes	Enabled
SpotDetection3D-IJ (1.7)	Last release	✓ Yes	Enabled
SpotDetection3D-Icy (v2.0)	Deprecated	✓ Yes	Disabled

Finish

- You can now use your workflow in BIAFLOWS
 - See how it performs (metrics)
 - Compare the results with other workflows
- Run the workflow locally on your machine, using docker to batch-analyze your images
- Link the executable workflow in your publication